



Read carefully before using the product. For use in “in vitro” diagnosis only

INTENDED USE

autoclean B is a solution used, once diluted, for internal washing of the hemostasis analyzer fluidic system.

COMPOSITION

Each autoclean B kit contains:

12 vials of 125 ml autoclean B concentrate (16x). The solution contains a tensioactive with colouring and sodium azide (0.1%, once diluted) as a preservative.

WARNING: The concentrated solution contains 1.6% sodium azide. Sodium azide may react with lead and copper plumbing to form highly explosive metal azide. On disposal, flush with a large volume of water to prevent azide buildup.

Preparation of washing solution



Dilute the contents of a 125 ml bottle of autoclean B (16x) with distilled or desmineralizada water to obtain a final volume of 2 litres.

STABILITY

The autoclean B concentrated solution is stable up to the expiry date stated on the label if stored at 18-25 °C.

Stability once diluted: At room temperature, it has the same stability as when unopened, provided that it is handled and stored correctly, avoiding contamination of the contents and prolonged exposure to temperatures that are not suitable for storage.

METHOD

autoclean B, once diluted, must be placed in the container of the hemostasis analyzer intended for this use (see the hemostasis analyzer user manual).

LIMITATIONS

Limitations in relation to the reagent:

Before use, inspect the reagent.

1. Do not use if the container is deteriorated and there is a loss of content.
2. Do not use if microbial or fungal contamination is observed.
3. Use before the expiry date.

WARNINGS AND PRECAUTIONS

autoclean B may suffer a loss of colour during the product life without affecting its functionality.

The use of “in vitro” diagnosis reagents for professional use requires taking into account the following indications:

- The product must only be used by qualified personnel.
- Once used, the product must be disposed of in special containers for biological waste.
- autoclean B (16x concentrate) contains sodium azide 1.6%.
R 23/24/25: Toxic by inhalation, in contact with skin and if swallowed.
R 32: Contact with acids liberates very toxic gas.
R 52/53: Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.
S 27: Take off immediately all contaminated clothing.
S 28: After contact with skin, wash immediately with plenty of water.
S 36/37: Wear suitable protective clothing and gloves.
S 45: In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
S 60: This material and its container must be disposed of as hazardous waste.
S 61: Avoid release to the environment. Refer to special instructions / safety data sheets.
- Diluted autoclean B contains sodium azide 0.1%.
R 20/21/22: Harmful by inhalation, in contact with skin and if swallowed.
S 36/37: Wear suitable protective clothing and gloves
- If you have doubts or need further information on the use of this product, consult the authorised distributor in your country.

PRESENTATION









215504 autoclean B 12 x 125 ml

Revision date: October, 2005

This document is available in several languages. The translations have been made from the master document in English. In the event of doubts or discrepancies, the wording in the master document in English shall take precedence.

Diagnostic Grifols, S.A. Passeig Fluvial, 24 - 08150 Parets del Vallès, ESPAÑA (SPAIN)



	T	Toxic
	N	Dangerous for the environment
	IVD	In vitro diagnostic medical device
	LOT	Batch code
		Use by
		Temperature limitation
		Consult Instructions for Use
	REF	Catalogue number